

	TECHNICAL SHEET	Document number: FT- BPV 02
	HIBANA®	Issue: Nov-11-2020
	<i>Korin Bio Inputs Line Activator</i>	Review: 01 Date: Dec-01-2020

ABOUT THE WORD HIBANA®

Within the Japanese tradition, fire is the beginning of everything.

It generates the earth starting the cycle of nature. The admiration for this element goes beyond, recalls the lantern festival that honors the dead and the Yamayaki - fire in the mountain that literally celebrates the end of millenary conflicts. Hibana means spark of life - is an allusion to this tradition and its power inducing the multiplication of life bringing productivity to your business.

BRAND

HIBANA® is an activator developed to facilitate and improve the activation of Korin bio-inputs.

COMPOSITION

Organic compost and minerals.

PRESENTATION

Physical nature: solid.

PACKAGES

15 kg bucket: black plastic packaging, with lid, capacity: 18 L, height 35,7 cm, diameter of mouth 30,1 cm, diameter of bottom 26,2 cm, net weight: 15 kg and gross weight: 15,7 kg.



Product specifications are on the label printed on the bucket.

INDICATION OF USE

HIBANA® is a single activator of Korin bio-inputs.

WAYS OF APPLICATION OF HIBANA®

Option A: application for activation of Korin bioinputs.

For activation of 100 L:

- 10 L of Korin bio-input;
- 1,5 kg of Hibana;
- 90 L of clean water, free of chlorine.

Instructions for preparation:

- Use only clean water, free of chlorine;
- Use hygienized containers, without chemical residues, specially bactericides and fungicides;
- Add the Hibana, already pre-diluted in the container with the clean water;

- Add Korin bio-input and homogenize;
- Cover the container so that no dust or insects enter;
- Rest for 48 to 72 hours in a cool and shady place and keep in them conditions until use;
- Do not agitate or aerate, activation should occur under low oxygenation;
- Use the activated solution in 15 days maximum.

Option B: direct application.

- Dilute the Hibana in a part of water and homogenize, directly in the spraying with Korin bio-input;
- For 10 L of Korin bio-input, use 1.5 kg of Hibana;
- Apply as recommended dosage of Korin bio-input.

STORAGE CONDITIONS

The storage must be at room temperature, in the shade, protected from direct sunlight.

DEACTIVATION

HIBANA[®] can be deactivated by any common sanitizing agent.

PRECAUTIONS

Avoid contact with eyes, skin and clothes.

To make the solution, do not use containers that were used for dilution of chemical products, such as defensives and disinfectants. We do not recommend the use of along with chemicals used as sanitizers (e.g. sodium hypochlorite).

PRODUCT EXEMPT FROM REGISTRATION

Distributed exclusively by:

KORIN AGRICULTURA E MEIO AMBIENTE LTDA.

Rodovia SP 191 | km 81,5
Estrada Municipal de Camaquã, s/nº
Bairro Área Rural | Ipeúna | SP | CEP.: 13537-000
CNPJ 08.183.347/0001-29
Tel.: (19) 35769518
www.korinagricultura.com.br
contato@korinagricultura.com.br
sac@korinagricultura.com.br

BRAZILIAN INDUSTRY

VALIDITY: 24 MONTHS